



OFFICE OF THE COLLECTOR AND DISTRICT MAGISTRATE, JAGATSINGHPUR
SOCIAL SECURITY SECTION

Tender Call Notice No 309

Dt: 04.02.2025

Sealed quotation/tenders are invited from the intending Registered firm/Supplying Agency for Supply of required items for smooth functioning of the District Disability Rehabilitation Centre (DDRC), Jagatsinghpur mentioned at Annexure-B for District Social Security Officer, Jagatsinghpur District.

Details of quotation/Tender papers along with the terms & conditions are available at this district website i.e jagatsinghpur.odisha.gov.in.

The quotation/tender paper (each page) must signed by the quotationer/tenderer as of acceptance of the term and Conditions.

The quotation/tender paper complete in all respect along with all required documents along with a Demand draft of Rs 20,000/- (Twenty Thousand) only from any nationalized bank towards cost of bid document in favour of the District Social Security Officer, Jagatsinghpur payable at Jagatsinghpur must be submitted in a sealed cover through Registered/Speed Post only to the District Social Security Officer, Jagatsinghpur, Collectorate, Jagatsinghpur, 754103 on or before 11 March, 2025 by 01.30 pm.

The quotation / tender shall be opened on 17 March, 2025 at 4.30 pm by the DMT (District Management Team) for DDRC in presence of the quotationer/tenderer or his/her authorised representatives only.

The quotation (s) / tender (s) received beyond the scheduled date and time will be automatically rejected.

The DMT reserves the right to accept/reject any or all quotation / tender paper without assigning any reason thereof.

Key/4/3/25
District Social Security Officer
Jagatsinghpur

Memo No 310 Dt 04.02.2025 //

Copy forwarded to the DIO, NIC, Jagatsinghpur for uploading the notice in the District official website of Jagatsinghpur for wide publicity.

Key/4/3/25
District Social Security Officer
Jagatsinghpur

Memo No 311 Dt 04.02.2025 //

Copy submitted to the ADM (General), Jagatsinghpur for favour of kind information & necessary action.

Key/4/3/25
District Social Security Officer
Jagatsinghpur


Memo No 312 Dt 04.02.2025 //

Copy submitted to PA to Collector, Jagatsinghpur for favour of kind information of the Collector & District Magistrate, Jagatsinghpur.

Key/4/3/25
District Social Security Officer
Jagatsinghpur

TERMS AND CONDITIONS FOR QUOTATIONERS/TENDERERS SUPPLY OF THE MATERIALS

1. The intending registered firm/Supplying agencies/Authorised Dealer must have valid GST,PAN and IT return certificate which Xerox copy should be attached with the quotation / tender paper.
2. The organisation should not have been black listed by the State/Central Govt or any other organisations.
3. Rate should be quoted inclusive of GST, Labour Cost for loading and unloading transportation charge and Installation charge of the materials.
4. The defective and damaged materials if found must be replaced by the agencies at their own cost within 3 days of report
5. The quotation / tender paper attached at Annexure-A should be filled in properly and eligibly without any correction/ overwriting.
6. The successful bidder will have to deliver the materials in the premises of the DDRC, Singharpur, ASRA School Campus,Nuapola,Jagatsinghpur. 754103 within 5 days from the date of issue of work order.
7. All materials will have to deliver to the above address on their own cost. The transposition charges, loading and unloading of materials, labour charges shall be borne by the supplier.
8. The bidder shall submit the quotation/tender documents within a sealed envelope super scribed as **“SUPPLY OF REQUIRED ITEMS FOR DDRC,JAGATSINGHPUR FOR DISTRICT SOCIAL SECURITY OFFICER,JAGATSINGHPUR”**.
9. The detail specification of the materials is attached at Annexure-B.
10. If any dispute arises the matter shall be referred to the DMT (District Management Team) of DDRC for settlement of the dispute.
11. The DMT reserves the right to reject any or all tenders without assigning any reason thereof and shall bear no liability whatsoever consequent upon decision.
12. The DMT will not responsible of missing or delay of any quotation/tender documents while in transit.


District Social Security Officer
Jagatsinghpur

Annexure-A

TENDER FORM

1. Name of the Bidder/Organisation :
(In Capitals Letter)
2. Complete Address of the Bidder/Organisation :
3. Name of the authorised signatory :
(In Capitals Letter)
4. Specimen Signature :
5. Mobile No :
6. GST Details :
7. Details of PAN/TIN or SRIN No :
(Attach self attested copy)
8. Total Bidding Price quoted :
9. Is the Bidder/Organisation is Black listed : Yes/No
By State/Central Govt(Attach Undertaking)
10. Details of IT return file :
(Attach Latest IT return File)

DECLARATION

I/We hereby certify that, the terms and conditions specification etc. Given with the tender notice has been read carefully and acceptable to me/us and that the information furnished above is complete and true to the best of my/our knowledge. I/We understand that in case of any deviation found in the above statement at any stage of the whole tender process then the tender paper submitted by me/us is liable to cancel. Further I/We certify that my/our Organisation is not have been black listed by any state/central Govt or any other organisation.

Place:

Date:

Signature with seal of the authorised signatory

24/4/25

Annexure-B

PRODUCT	SPECIFICATION
SUPINATOR PRONATION MACHINE	<p>It is used for supination and pronation exercises of the wrist. Two precision chrome plated circular cylindrical frame with a grip handle are mounted on a laminated wooden board. The resistance of these frames are adjusted with help of rotating knob, which adjusts the tensions of internal spring assembly. It is used for supination and pronation exercises of the wrist. Two precision chrome plated circular cylindrical frame with a grip handle are mounted on a laminated wooden board. The resistance of these frames are adjusted with help of rotating knob, which adjusts the tensions of internal spring assembly.</p>
PEANUT BALL	<p>Exercise with peanut yoga ball can relax your muscle, relief back pain cultivate a peaceful center and a sense of balance. It is great for fitness, therapy health, Labor Birthing exercise, kid's play and dog training. Constructed from durable and thick resilient PVC material that makes it Anti-Burst tested and Non-Slip. The shape of a peanut allow it only moves in one direction offering a stable base for balance. Includes a quick inflation foot pump. It take a few minutes to completely inflate a ball. It can use at the gym, in yoga, or in your living room. peanut ball size 45*90cm</p>
STATIC CYCLE	<p>Cardio Benefits Riding stationary bicycle counts as cardiovascular exercise. Regular cardio exercise helps lower high blood pressure, regulates blood sugar and may prevent heart attacks, especially if you've already suffered from one. Regularly riding a stationary bike helps your heart become more efficient at pumping blood. It also increases the amount of good cholesterol while lowering the amount of bad cholesterol in your system to result in less plaque build-up in your arteries. The simple act of riding a stationary bicycle 30 minutes a day, five times per week along with proper diet will help you to actively manage a healthy lifestyle. Joint-Friendly Many types of aerobic activity, such as jogging or team sports, can be hard on your joints because of the impact involved. A stationary cycle enables you to get your heart pumping without putting undue stress on your joints. The bike is especially good for people as it reduces mechanical stress on the back, hips, knees and ankles even when compared to walking. You can indulge in daily exercise with a lower risk of injury. Installation, Warranty and After Sales Support : To help our customers in product installation we provide an installation manual in the product box. YouTube installation video is available on-demand and an introduction to third-party technicians who provide installation on chargeable basis. In case there is any damage doing transit then you can write to us sharing the images of the defective product damaged parts will be replaced free of cost</p>
BALANCE BOARD	<p>Balance Boards are being used in physical education classes around the globe to help children to develop not only their balance and coordination, but also their posture and core strength. They are also proving effective to help children with conditions such as cerebral palsy to improve their coordination.</p>

Reshmi/25

INFRARED LAMP	Infrared lamps are extensively used for heating purposes. Our infrared lamps are both domestically and commercially used owing to its superior durability and excellent performance. Provides excellent heating and illumination Compact design Easy to install and use Long lasting service life Cost efficient
MYOFASCIAL RELEASE TRIGGER POINT THERAPY TOOL	
MIRROR THERAPY	Mirrors Therapy are used for training individuals in walking after a stroke or other affections that force them to relearn or correct their walking or posture. Mirrors provide Visual feedback of patients progress.
ISTAM TOOL	Refers to a technique that uses instruments to remove scar tissues from injured soft tissues and facilitate healing process through formation of new extracellular matrix proteins such as collagen. Recently, frequent use of this instrument has increased in the fields of sports rehabilitation and athlete training.
TAPING	
SHORTWAVE DIATHERMY	500 is an advance device, unique design, high power shortwave frequency .The high frequency generating circuitry produces radio frequency waves of 27.12 MHz that penetrates deep into muscle tissue, with negligible heating in fatty and hard tissue. The unit have facility to produce the controlled energy into the area to be treated with marginal effect on adjacent tissues..Model HME-001 Operating Frequency 27.12 MHz Wave Length 11.02 meter RF Power Source Vacuum Generating Output power 500 watts Mode of operation Continuous and pulse Pulse Frequency 2-200 Hz Variable Fuse 10 Amps Power Display Analogue display Treatment timer 0-99 minutes. Weight 38 Kg
CUPPING KIT	Kinesio Tex Tape is used to improve joint alignment by affecting the muscles and fascia and can reduce poor function of a joint by influencing opposing muscle groups and joint mobility.
RAMP WALKING	
CPM MACHINE	Knee Continuous Passive Motion unit cpm knee Physical Therapy is a unique device to improve the flexion movement; the stiff muscle do not allow the full range of motion due to less movement, exercise of joints due to pain after surgical precedure, accidental injury Polio, paralysis and hemiplegia etc. - 132 is a micro computer based, technology derived dielectro mechanical device which is safe and reliable to use. It constantly moves the joint through a controlled range of motion; the exact range is dependent upon the joint. Rehabilitation : control post-operative pain, reduce inflammation, provide passive motion in a specific plane of movement, and protect the healing repair or tissue. Technical Specifications : Input Voltage : 220v; Treatment timer : 0-99 min, Power Consumption : 35 w max., Flexions hold time : 0-9 sec.; Extension hold time : 0-9sec.; Flexsion Angle : 5 - 120 degree.

Regy. 3.25

<p>PHYSIOTHERAPY EXERCISE CHAIR</p>	<p>Effective quadriceps and hamstrings exerciser without changing the position of patient. The unit consists of foam padded table and a torque unit with two lever arms, one of which holds the weight while the other provides contact with the patient. By changing the angle between the two arms, maximum or minimum resistance can be given at any point in the range. Complete with six weight plates of 1 kg each.</p>
<p>SHOCK WAVE THERAPY</p>	<ul style="list-style-type: none"> • Can be used to treat a wide variety of musculoskeletal conditions, in particular, those involving areas where major connective tissue attaches to bone. Common sites that can be successfully treated are: <ul style="list-style-type: none"> • Foot - plantar fasciitis and Achilles tendonitis • Knee - patellar tendonitis and shin pain • Elbow - tennis or golfer's elbow • Shoulder - rotator cuff tendonitis and calcific • Hip - trochanteric bursitis • Muscles - various trigger points throughout body • EDSWT (Erectile Dysfunction Shock Wave Therapy) is a new noninvasive therapy approach using low intensity Extracorporeal Shockwave technology for the treatment of Erectile Dysfunction. • Energy : 10mj (0.2Bar) to 190mj (5Bar). adjustable stepping 10mj • Frequency : 1HZ to 22 HZ adjustable stepping at 1HZ • Transmitter : 5pcs • Operating mode : Continuous/4/8/16 pulses • Input Power:- 220 -240 V ~; 50/60 Hz <p>Power:140W max</p>

20/4/3/25

JSBHF 14 WHOLE BODY EXERCISER	<p>Technical specification</p> <p>Innovative body shaping equipment that effectively firms and tones your body. Designed with resistance cord system for more intense workout. Enhance workout intensity by widening stepping angels to stimulate more muscle movement. Auto programs (circulation, Firm and Strengthen) that integrates different speeds for optimal workout. Manual program with adjustable speed and timer to suit individual workout preference. Wheel caster for easy maneuverability. Adjustable rubber sands help to stabilize the device on an uneven surface. Designed with motor protecting system to prevent event of improper operation such as overloading and uneven weight distribution at one side of the vibrating platform.</p>
4D ELECTRIC MASSAGE MAT	<p>The "4D" of these massage chairs refers to advanced massage roller capability. Users get the benefit of having massage rollers or nodes that mimic the action of human hands, allowing for massaging in directions that standard massage chairs simply can't manage themselves.</p>
ELECTROSURGICAL CAUTERY MACHINE	
ELECTRICAL STIMULATOR	<p>Technical Specifications Input Voltage : 220 V AC 50 Hz Galvanic : Maximum output 100 Volt. Interrupted galvanic : A switch to select 4 different mode of Frequency. Faradic : Compriser 40 pulse of 0.9 ms. each pulse of 0.9 ms. Surge Faradic : Step less control to increases.</p>

02/04/2025

**INTERFERENTIAL THERAPY
(IFT)**

INTERFERENTIAL THERPY 70 program

TECHNICAL SPECIFICATIONS

- ▶ Therapy mode: 4P, 4 PV, 2P, Russian current
- ▶ Vector scan 4 P 90°, 2P 45°
- ▶ 70 pre-program with protocols
- ▶ Spectrum pattern: Triangular A, Trapezoidal Square
- ▶ Treatment Duration:-1-99 min in 1 min increments
- ▶ Alarms Audio Alarm at the end of Treatment with auto cut off
- ▶ Beat frq: 0-150cont, 0-10 scan, 80-150 scan & 150 scan
- ▶ Carrier frq: 2 khz & 4 khz, 2.5 khz
- ▶ Micro computer controlled based with baeklight 16×2 alpha numeric LCD Display of timer, therapy mode, carrier frequency makes, Separate display for both channel out put amplitude, Beat frq
- ▶ Amplitude controlled :- 0 to 99 ma (peak to peak variable)
- ▶ Amplitude balance controlled: To equilibise 2 channel current
- ▶ Input Power:- 220-240 V ~, 50/60 Hz
- ▶ Weight: 4 kg approx
- ▶ Brand: physioift

STANDARD ACCESSORIES

- ▶ Main Unit: 01 No.
- ▶ Ift carbon electrode: 2 red & 2 blue
- ▶ Electrode sponge cover: 4 No
- ▶ Ift out put lead: 2 red wire & 2 black wire
- ▶ Fixing elastic belts: 02 set. (each set three different size of belt small, medium, & large)
- ▶ Power cable: 01 No.
- ▶ Instruction Manual :01 No
- ▶ Carry bag: 1 No
- ▶ .Jel bottiale: INo

copy 3.25

COLD THERAPY MACHINE	<p>SPECIFICATIONS:</p> <ul style="list-style-type: none"> • Construction: Unit's outer body made of 20 gauge and inner chamber of 18 gauges S.S. Sheet, double walled unit With well insulated Top & Sides. Covering with air tight painted steel sheet with high quality • Dimensions: Inner Tank size 40 x 30 x 35cm deep. : Over all, approx.60cm x 60cm x 102cm. • Temp Control: Fitted with Digital Temperature Indicator cum Controller to control Temperature to -12 to -6 deg Cen. : Fitted with 2 indicator Lamps. • Cooling Unit: Heavy Duty Sealed Unit of "Kirloskar" or "Tecumseh" makes. • Cold Packs. : 6 Packs of Standard size (28x36cm) • Drain Valve: A drain valve is provided for easy defrosting and cleaning. • Mobility: Mounted on Four, 75mm ball bearing castors for easy mobility. • Accessories: Six S.S Trays are provided to keep Cold packs in the unit. • Finish: Polished finish inside and special paint finish outside. Power: 220-240 V AC.
CONTRAST BATH THERAPY	<p>Technical Specification</p> <ul style="list-style-type: none"> • Made of well insulated water tanks of 20 gauges S.S 304 Grade outer chamber CRC with power coating paint & Covering with painted steel sheet. Tanks are covered with laminated wooden rim & covers. • Energy saving glass wool insulation • Water Tanks: Two, Size: 50cm x 38cm x 40cm deep. • Hot-Bath Tank: Fitted with One, 1kw special heating element. • Cold-Bath Tank: Fitted with heavy duty cooling compressor unit • Temperature Control: Each tank is provided with Digital Thermometer, Thermostat for Temperature Control. & Pilot lights. • Mounting: Mounted on 10cm dia. Four casters. • Drainage: Drainage system for emptying. Power: 220-240 V AC.

Ref 4.3.25

<p>CERVICAL CUM LUMBER TRACTION MACHINE</p>	<p>COMPUTERIZED DIGITAL TRACTION MACHINE (soft touch) that are widely used across the medical industry for delivering accurate results. These instruments are made from authentic material so as these instruments are safe to use and are based on contemporary technology. Our instruments are a classic example of the technological advancements that have taken place in the medical sector</p> <p>TRACTION TABLE: Height 2.5 Feet Weight Upto 70 Kg Body Material SS Usage Clinical, Hospital</p>
<p>HYDROCOLLATOR</p>	<p>Hydrocollator is used for heating the steam pack .Technical Specification: Inner Chamber Size:17" X 13" X (14"H) Outer Chamber Size: 20" X 16" X (19"H) Heater, 4 PACK.</p>
<p>SHOULDER WHEEL EXERCISER MAGNETIC</p>	<p>Shoulder wheel exerciser magnetic</p> <p>TECHNICAL SPECIFICATIONS</p> <ul style="list-style-type: none"> ▶ 24" Inch Handle Rod for Smooth Shoulder Rotation, & for Light weight, easy to install ▶ Fully automatic 5-function digital LCD display: speed, elapsed time, calories, distance and scan of all functions ▶ Quiet and smooth adjustable Resistance is controlled by turning the resistance knob. ▶ shoulder wheel height adjustable vertical wallmount direction

copy/4/3/25